


Safety Communication

رسالة سلامة

Y-piece and the associated adhesive joint did not conform to specification

Device/ Product Description:	Breathing Circuits, Ventilator
Affected Products:	<ul style="list-style-type: none"> • MP00333 VentStar Anesthesia (N) 180 • MP00353 VentStar Basic (N) 180 • MP00363 VentStar Watertrap (N) 180 • MP01327 ID Circuit Basic (N) 180 • MP01328 ID Circuit Watertrap (N) 180
Manufacturer:	Draeger Medical Systems Inc
Problem:	The Y-pieces, due to a manufacturing error, can in some cases become disconnected at tensile forces below what is specified. If the Y-piece becomes detached from the hose during use, a leak will be caused that can hinder ventilation of the patient.
Recommendation /Actions:	<ul style="list-style-type: none"> • Make sure that this document, with accompanying notice is reached to the end-users. • Check the quality of the adhesive joint before use by holding the Y-piece with one hand and then pulling (not twisting) each of the two hoses one after another so that the hose just starts to stretch visibly. • Use the hoses with devices that manufacturer has declared compatible in the Instructions for Use or accessories list of the basic device. • Contact the authorized representative for required support, and to replace the unused stocks that were manufactured between 01.01.2018 and 13.11. 2020.. <p>For more information, Please click here.</p> <p>If you think you had a problem with your device or a device your patient uses, please report the problem to SFDA through: NCMDR Vigilance system</p>

	(19999) unified call center	
Devices/Products photo:		
Authorized Representative Details	AR name:	Draeger Arabia Co Ltd
	Assigned Contact Person:	Abdulrahman Alsaeed
	Mobile/Phone:	+(966) 555556548
	Email:	abdulrahman.alsaeed@draeger.com