


## Safety Communication

## رسالة سلامة

### Potential risk when using the system with specific software versions

<b>Device/ Product Description:</b>	In-vitro diagnostics - microbiological products	
<b>Affected product:</b>	Model Name:	Accelerate Pheno system.
	Catalogue No.:	10401008, 10301008
<b>Manufacturer:</b>	Accelerate Diagnostics Inc	
<b>Problem:</b>	Potential risk of very major discrepancies, major discrepancies and/or essential agreement <90% with some organism-antimicrobial combinations as listed below, when using the Accelerate Pheno™ system with software versions 1.4.1 and later.	
<b>Recommendation /Actions:</b>	<ul style="list-style-type: none"> <li>• Make sure that this document is reached to the end-users.</li> <li>• Follow the instructions as indicated in the attached Field Safety Corrective Action.</li> <li>• Contact the authorized representative for required corrective action.</li> </ul> <p>For more information please click <a href="#">here</a>.</p> <p>If you think you had a problem with your device or a device your patient uses, please report the problem to SFDA through:  <a href="#">NCMDR</a>  <a href="#">Vigilance system</a>            (19999 )unified call center</p>	

<b>Devices/Products photo:</b>		
<b>Authorized Representative Details</b>	<b>AR name:</b>	Creative Healthcare Establishment
	<b>Assigned Contact Person:</b>	Khaled Alressini
	<b>Mobile/Phone:</b>	+(966) 503478751
	<b>Email:</b>	info@chc.med.sa

